





RAF1 Antibody

Product Code	CSB-PA080220
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04049
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC; WB:1:1000-2000, IHC:1:50-200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	RAF proto-oncogene serine/threonine-protein kinase
Immunogen Species	Homo sapiens (Human)
Target Names	RAF1

Image



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma using Raf-1 Polyclonal Antibody.



Western blot analysis of Hela using Raf-1 Polyclonal Antibody.