





Pan Methylated Lysine Monoclonal Antibody

Product Code	CSB-MA080224
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Conjugated Protein
Raised In	Mouse
Species Reactivity	All
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:50-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Product Type	Monoclonal Antibody
Immunogen Species	All
Clone No.	7H8

Image



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using Pan Methylated Lysine Monoclonal Antibody.



Western blot analysis of Hela using Pan Methylated Lysine Monoclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.