





MAP2K5 Antibody

Product Code	CSB-PA080231
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13163
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC; WB:1:1000-2000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Dual specificity mitogen-activated protein kinase kinase 5
Immunogen Species	Homo sapiens (Human)
Target Names	MAP2K5

Image



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using MEK-5 Polyclonal Antibody.



Western blot analysis of 1) HeLa, 2) 3T3, 3) PC12 using MEK-5 Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.