




HAO1 Antibody

Product Code	CSB-PA080232
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJM8
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Mouse,Rat
Tested Applications	WB,IHC;WB:1:2000-5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Hydroxyacid oxidase 1
Immunogen Species	Homo sapiens (Human)
Target Names	HAO1

Image



Western blot analysis of 1) Mouse Liver Tissue, 2) Rat Liver Tissue using HAO1 Polyclonal Antibody.