




HAO1 Antibody

Product Code	CSB-PA080232
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJM8
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Mouse,Rat
Tested Applications	WB,IHC;WB:1:2000-5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Hydroxyacid oxidase 1
Immunogen Species	Homo sapiens (Human)
Target Names	HAO1

Image



Western blot analysis of 1) Mouse Liver Tissue, 2) Rat Liver Tissue using HAO1 Polyclonal Antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.