





TGFB1 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA080233 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P01137 |
| Immunogen | Recombinant Protein |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | WB, IHC; WB:1:1000-2000, IHC:1:200-500 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Isotype | IgG |
| Alias | Transforming growth factor beta-1 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | TGFB1 |

Image



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using TGFB1 Polyclonal Antibody.



Western blot analysis of 1) Mouse Liver Tissue, 2) Rat Liver Tissue using TGFB1 Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.