





TGFB1 Antibody

Product Code	CSB-PA080233
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01137
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC; WB:1:1000-2000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Transforming growth factor beta-1
Immunogen Species	Homo sapiens (Human)
Target Names	TGFB1

Image



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using TGFB1 Polyclonal Antibody.



Western blot analysis of 1) Mouse Liver Tissue, 2) Rat Liver Tissue using TGFB1 Polyclonal Antibody.