





STAT1 Antibody

Product Code	CSB-PA080236
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42224
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC; WB:1:1000-2000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Signal transducer and activator of transcription 1-alpha/beta
Immunogen Species	Homo sapiens (Human)
Target Names	STAT1

Image



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using Stat1 Polyclonal Antibody.



Western blot analysis of 1) HeLa, 2) 3T3, 3) Rat Heart Tissue using Stat1 Polyclonal Antibody.