




BSA Monoclonal Antibody

Product Code	CSB-MA080239
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Purified Protein
Raised In	Mouse
Species Reactivity	Cow
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Product Type	Monoclonal Antibody
Immunogen Species	Cow
Clone No.	9E7

Image



Western blot analysis of Bovine Serum Albumin using Bovine Serum Albumin Monoclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.