





GAPDH Antibody

Product Code	CSB-PA080240
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04406
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC; WB:1:5000-10000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Glyceraldehyde-3-phosphate dehydrogenase
Immunogen Species	Homo sapiens (Human)
Target Names	GAPDH

Image



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using GAPDH Polyclonal Antibody.



Western blot analysis of HeLa using GAPDH Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.