





ACTN1 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA080241 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P12814 |
| Immunogen | Recombinant Protein |
| Raised In | Rabbit |
| Species Reactivity | Mouse,Rat |
| Tested Applications | WB, IHC; WB:1:1000-2000, IHC:1:200-500 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Isotype | IgG |
| Alias | Alpha-actinin-1 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | ACTN1 |

Image



Immunohistochemical analysis of paraffin-embedded Mouse Kidney Tissue using Actinin- α 1 Polyclonal Antibody.



Western blot analysis of 1) 3T3, 2) PC12, 3) Rat Heart Tissue using Actinin- α 1 Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.