





# STAT5A Antibody

<b>Product Code</b>	CSB-PA080243
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P42229
<b>Immunogen</b>	Recombinant Protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	WB, IHC; WB:1:500-1000, IHC:1:100-200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Isotype</b>	IgG
<b>Alias</b>	Signal transducer and activator of transcription 5A
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	STAT5A

## Image



Immunohistochemical analysis of paraffin-embedded Human Ovarian Carcinoma using Stat5a Polyclonal Antibody.



Western blot analysis of K562 using Stat5a Polyclonal Antibody.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.