





AKT2 Antibody

Product Code	CSB-PA080247
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31751
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC; WB:1:500-1000, IHC:1:100-200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	RAC-beta serine/threonine-protein kinase
Immunogen Species	Homo sapiens (Human)
Target Names	AKT2

Image



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma using Akt2 Polyclonal Antibody.



Western blot analysis of 1) HeLa, 2) 3T3, 3) PC12 using Akt2 Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.