





CDKN1A Antibody

Product Code	CSB-PA080248
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38936
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC; WB:1:1000-2000, IHC:1:100-200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	Cyclin-dependent kinase inhibitor 1
Immunogen Species	Homo sapiens (Human)
Target Names	CDKN1A

Image



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using p21 Polyclonal Antibody.



Western blot analysis of 1) 293T, 2) MCF7, 3) Mouse Liver Tissue, 4) Rat Liver Tissue using p21 Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.