




HSPB8 Antibody

Product Code	CSB-PA080249
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJY1
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Heat shock protein beta-8
Immunogen Species	Homo sapiens (Human)
Target Names	HSPB8

Image



Western blot analysis of Hela using Crystallin- α C Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.