



SMAD2 Antibody

Product Code	CSB-PA080253
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15796
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB;WB:1:1000-2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	Mothers against decapentaplegic homolog 2
Immunogen Species	Homo sapiens (Human)
Target Names	SMAD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.