



TUBB Monoclonal Antibody

Product Code	CSB-MA080255
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07437
Immunogen	Recombinant Protein
Raised In	Mouse
Species Reactivity	H,R,M,Mk,Dg,C,Hm,Rb,Sh,Insect,Yeast
Tested Applications	ELISA,WB,IHC;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500
Form	Liquid
Conjugate	HRP
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	β-tubulin(HRP Conjugated)
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Gene Names	TUBB