



HSPB1 Monoclonal Antibody

Product Code	CSB-MA080263
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04792
Immunogen	Recombinant Protein
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB: 1:500-1:5000, IHC: 1:50-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	HSP27
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HSPB1
Clone No.	16G7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.