





Phosphotyrosine Monoclonal Antibody

Product Code	CSB-MA080265
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	N/A
Immunogen	Purified Protein
Raised In	Mouse
Species Reactivity	All
Tested Applications	ELISA,WB,IHC;Recommended dilution:WB:1:500-1:2000,IHC:1:100-1:300
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	Phosphotyrosine
Product Type	Monoclonal Antibody
Clone No.	13B7

Image



Western blot analysis of 1) Jurkat, 2) 293T, 3) Rat Brain Tissue, 4) Mouse Brain Tissue with Phosphotyrosine Mouse mAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using Phosphotyrosine Mouse mAb diluted at 1:200.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.