




TNF Antibody

Product Code	CSB-PA080267
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01375
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC; WB:1:500-1:2000, IHC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	TNF α
Immunogen Species	Homo sapiens (Human)
Target Names	TNF

Image



Western blot analysis of SH-SY5Y 293T 3T3 lysis using TNF α antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.