



IGFN1 Antibody

Product Code	CSB-PA080268
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86VF2
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Pg
Tested Applications	WB, IHC; WB:1:500-1:2000, IHC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	eEF1A2 binding protein
Immunogen Species	Homo sapiens (Human)
Target Names	IGFN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.