



TNF Monoclonal Antibody

Product Code	CSB-MA080272
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01375
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	TNF α
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TNF
Clone No.	9F10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.