



Mono-Methyl-TP53(Lys370) Monoclonal Antibody

Product Code	CSB-MA080273
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04637
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC;Recommended dilution:IHC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	P53(Mono Methyl Lys370)
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TP53
Clone No.	25G1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.