



FLOT2 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA080278 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q14254 |
| Immunogen | Recombinant Protein |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | WB;WB:1:500-1:2000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Isotype | IgG |
| Alias | Flotillin-2 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | FLOT2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |