



BCL2L1 Antibody

Product Code	CSB-PA080279
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07817
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	IHC;IHC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	Bcl-XI
Immunogen Species	Homo sapiens (Human)
Target Names	BCL2L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.