



BEST1 Antibody

Product Code	CSB-PA080283
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O76090
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB;WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	Bestrophin-1
Immunogen Species	Homo sapiens (Human)
Target Names	BEST1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.