



CACNA2D1 Antibody

Product Code	CSB-PA080288
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54289
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	WB;WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	CaV α 2 δ 1
Immunogen Species	Homo sapiens (Human)
Target Names	CACNA2D1