



HSP90AA1 Monoclonal Antibody

Product Code	CSB-MA080290
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07900
Immunogen	Recombinant Protein
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	HSP90 α
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Gene Names	HSP90AA1