



HSPA8 Monoclonal Antibody

| | |
|----------------------------|--|
| Product Code | CSB-MA080291 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P11142 |
| Immunogen | Recombinant Protein |
| Raised In | Mouse |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB;Recommended dilution:WB:1:500-1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol. |
| Isotype | IgG |
| Alias | HSC70 |
| Product Type | Monoclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | HSPA8 |
| Clone No. | 17B7 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |