




H2AFZ Antibody

Product Code	CSB-PA011027
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C0S5
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Mouse,Rat
Tested Applications	WB;WB:1:1000-2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	H2AFZ; H2AZ; Histone H2A.Z; H2A/z
Immunogen Species	Homo sapiens (Human)
Target Names	H2AFZ

Image



Western blot analysis of 1) Rat Brain, 2) Mouse Brain tissue, diluted at 1:2000.