




SNRPN Antibody

Product Code	CSB-PA011056
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P63162
Immunogen	Synthesized peptide derived from the Internal region of Human SNRPN.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SNRPN; HCERN3; SMN; Small nuclear ribonucleoprotein-associated protein N; snRNP-N; Sm protein D; Sm-D; Sm protein N; Sm-N; SmN; Tissue-specific-splicing protein
Immunogen Species	Homo sapiens (Human)
Target Names	SNRPN

Image



Western Blot analysis of HeLa cells using SNRPN Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.