




IL15RA Antibody

Product Code	CSB-PA011216
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13261
Immunogen	Synthesized peptide derived from the Internal region of Human IL-15R α .
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	IL15RA; Interleukin-15 receptor subunit alpha; IL-15 receptor subunit alpha; IL-15R-alpha; IL-15RA; CD antigen CD215
Immunogen Species	Homo sapiens (Human)
Target Names	IL15RA

Image



Western Blot analysis of COLO cells using IL-15R α Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.