



Phospho-AR (S213) Antibody

Product Code	CSB-PA011251
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10275
Immunogen	Synthesized peptide derived from Human AR around the phosphorylation site of S213.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC,ELISA;IHC:1/100-1/300.ELISA:1/40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	AR; DHTR; NR3C4; Androgen receptor; Dihydrotestosterone receptor; Nuclear receptor subfamily 3 group C member 4
Species	Homo sapiens (Human)
Target Names	AR