



Phospho-AR (S650) Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA011252 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P10275 |
| Immunogen | Synthesized peptide derived from Human AR around the phosphorylation site of S650. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | IHC, ELISA; IHC:1:100-1:300, ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | AR; DHTR; NR3C4; Androgen receptor; Dihydrotestosterone receptor; Nuclear receptor subfamily 3 group C member 4 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | AR |