




AR Antibody

Product Code	CSB-PA011255
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10275
Immunogen	Synthesized peptide derived from Human AR around the non-phosphorylation site of S650.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	AR; DHTR; NR3C4; Androgen receptor; Dihydrotestosterone receptor; Nuclear receptor subfamily 3 group C member 4
Immunogen Species	Homo sapiens (Human)
Target Names	AR

Image



Western Blot analysis of 823 cells using AR Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.