



AR Antibody

Product Code	CSB-PA011257
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10275
Immunogen	Synthesized peptide derived from Human AR around the non-phosphorylation site of S94.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	AR; DHTR; NR3C4; Androgen receptor; Dihydrotestosterone receptor; Nuclear receptor subfamily 3 group C member 4
Immunogen Species	Homo sapiens (Human)
Target Names	AR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.