




AR Antibody

Product Code	CSB-PA011258
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10275
Immunogen	Synthesized peptide derived from Human AR around the non-phosphorylation site of Y363.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	AR; DHTR; NR3C4; Androgen receptor; Dihydrotestosterone receptor; Nuclear receptor subfamily 3 group C member 4
Immunogen Species	Homo sapiens (Human)
Target Names	AR

Image



Western Blot analysis of 293T NIH-3T3 cells using AR Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.