



Phospho-ITGA4 (S1027) Antibody

Product Code	CSB-PA011259
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13612
Immunogen	Synthesized peptide derived from Human Integrin α 4 around the phosphorylation site of S1027.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ITGA4; CD49D; Integrin alpha-4; CD49 antigen-like family member D; Integrin alpha-IV; VLA-4 subunit alpha; CD antigen CD49d
Immunogen Species	Homo sapiens (Human)
Target Names	ITGA4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.