



ITPKB Antibody

Product Code	CSB-PA011279
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27987
Immunogen	Synthesized peptide derived from the Internal region of Human InsP 3-kinase B.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ITPKB; Inositol-trisphosphate 3-kinase B; Inositol 1; 4,5-trisphosphate 3-kinase B; IP3 3-kinase B; IP3K B; InsP 3-kinase B
Immunogen Species	Homo sapiens (Human)
Target Names	ITPKB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.