




JAK3 Antibody

Product Code	CSB-PA011293
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52333
Immunogen	Synthesized peptide derived from Human JAK3 around the non-phosphorylation site of Y785.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	JAK3; Tyrosine-protein kinase JAK3; Janus kinase 3; JAK-3; Leukocyte janus kinase; L-JAK
Immunogen Species	Homo sapiens (Human)
Target Names	JAK3

Image



Western Blot analysis of 293T cells using JAK3 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.