



NDUFS7 Antibody

Product Code	CSB-PA011306
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75251
Immunogen	Synthesized peptide derived from the C-terminal region of Human NDUFS7.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NDUFS7; NADH dehydrogenase [ubiquinone] iron-sulfur protein 7; mitochondrial; Complex I-20kD; CI-20kD; NADH-ubiquinone oxidoreductase 20 kDa subunit; PSST subunit
Immunogen Species	Homo sapiens (Human)
Target Names	NDUFS7