



SFT2D2 Antibody

Product Code	CSB-PA010504
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95562
Immunogen	Synthesized peptide derived from the Internal region of Human SFT2D2.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SFT2D2; Vesicle transport protein SFT2B; SFT2 domain-containing protein 2
Immunogen Species	Homo sapiens (Human)
Target Names	SFT2D2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.