



Phospho-KCND2 (S616) Antibody

Product Code	CSB-PA010505
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZV8
Immunogen	Synthesized peptide derived from Human Kv4.2 around the phosphorylation site of S616.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	KCND2; KIAA1044; Potassium voltage-gated channel subfamily D member 2; Voltage-gated potassium channel subunit Kv4.2
Immunogen Species	Homo sapiens (Human)
Target Names	KCND2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.