




PTPRC Antibody

Product Code	CSB-PA010546
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08575
Immunogen	Synthesized peptide derived from Human CD45 around the non-phosphorylation site of S1007.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PTPRC; CD45; Receptor-type tyrosine-protein phosphatase C; Leukocyte common antigen; L-CA; T200; CD antigen CD45
Immunogen Species	Homo sapiens (Human)
Target Names	PTPRC

Image



Western Blot analysis of Jurkat cells using CD45 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.