



Phospho-RELB (S552) Antibody

Product Code	CSB-PA010548
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01201
Immunogen	Synthesized peptide derived from Human RelB around the phosphorylation site of S552.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RELB; Transcription factor RelB; I-Rel
Immunogen Species	Homo sapiens (Human)
Target Names	RELB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.