





RAD17 Antibody

Product Code	CSB-PA019253EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75943
Immunogen	Recombinant Human Cell cycle checkpoint protein RAD17 protein (1-270AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell cycle checkpoint protein RAD17 (hRad17) (RF-C/activator 1 homolog), RAD17, R24L
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RAD17

Image



Western blot
All lanes: RAD17 antibody at 12µg/ml + Mouse thymus tissue
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 78, 76, 58, 67 kDa
Observed band size: 78 kDa



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA019253EA01HU at dilution of 1:100