



BMPR1A Antibody, FITC conjugated

Product Code	CSB-PA002748LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36894
Immunogen	Recombinant Human Bone morphogenetic protein receptor type-1A protein (177-532AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bone morphogenetic protein receptor type-1A (BMP type-1A receptor) (BMPR-1A) (EC 2.7.11.30) (Activin receptor-like kinase 3) (ALK-3) (Serine/threonine-protein kinase receptor R5) (SKR5) (CD antigen CD292), BMPR1A, ACVRLK3 ALK3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	BMPR1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.