



RNF133 Antibody, Biotin conjugated

Product Code	CSB-PA855485LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVZ7
Immunogen	Recombinant Human E3 ubiquitin-protein ligase RNF133 protein (65-167AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase RNF133 (EC 2.3.2.27) (RING finger protein 133) (RING-type E3 ubiquitin transferase RNF133), RNF133
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RNF133
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.