



PPME1 Antibody, Biotin conjugated

Product Code	CSB-PA018501LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y570
Immunogen	Recombinant Human Protein phosphatase methylesterase 1 protein (201-386AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein phosphatase methylesterase 1 (PME-1) (EC 3.1.1.89), PPME1, PME1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPME1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.