





STXBP6 Antibody

Product Code	CSB-PA850881LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFX7
Immunogen	Recombinant Human Syntaxin-binding protein 6 protein (2-210AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Syntaxin-binding protein 6 (Amisyn), STXBP6
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	STXBP6


Image




Western blot
All lanes: STXBP6 antibody at 2μg/ml + Mouse brain tissue
Secondary
Goat polyclonal to rabbit at 1/10000 dilution
Predicted band size: 24, 12 kDa
Observed band size: 24 kDa



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA850881LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA850881LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA850881LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)